

## Myosin light chain 1/3, skeletal muscle isoform

Catalog No: #AP85899



Package Size: #AP85899-1 50ug #AP85899-2 100ug #AP85899-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Myosin light chain 1/3, skeletal muscle isoform
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MLC1/MLC3,Myosin light chain alkali 1/2,Myosin light chain A1/A2,Myf
Accession No.	Uniprot ID:P05977
Uniprot	P05977
GeneID	17901;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Myf1

Note: This product is for in vitro research use only