

Nuclear autoantigenic sperm protein

Catalog No: #AP85900



Package Size: #AP85900-1 50ug #AP85900-2 100ug #AP85900-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nuclear autoantigenic sperm protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	NASP
Accession No.	Uniprot ID:Q2T9P4
Uniprot	Q2T9P4
GeneID	512427;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NASP

Note: This product is for in vitro research use only