

Neuromedin-K receptor

Catalog No: #AP85908



Package Size: #AP85908-1 50ug #AP85908-2 100ug #AP85908-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neuromedin-K receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	NKR,NK-3 receptor,NK-3R,Neurokinin B receptor,Tachykinin receptor 3,Tac3r
Accession No.	Uniprot ID:P47937
Uniprot	P47937
GeneID	21338;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tacr3

Note: This product is for in vitro research use only