

## Oxysterols receptor LXR-beta

Catalog No: #AP85909



Package Size: #AP85909-1 50ug #AP85909-2 100ug #AP85909-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Oxysterols receptor LXR-beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Liver X receptor beta,Nuclear receptor NER,Nuclear receptor subfamily 1 group H member 2,Ubiquitously-expressed nuclear receptor,LXRB,NER,UNR
Accession No.	Uniprot ID:P55055
Uniprot	P55055
GeneID	7376;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:NR1H2

Note: This product is for in vitro research use only