

Nuclear receptor subfamily 1 group I member 2

Catalog No: #AP85910



Package Size: #AP85910-1 50ug #AP85910-2 100ug #AP85910-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nuclear receptor subfamily 1 group I member 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Orphan nuclear receptor PXR,Pregnane X receptor,Pxr
Accession No.	Uniprot ID:Q9R1A7
Uniprot	Q9R1A7
GeneID	84385;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nr1i2

Note: This product is for in vitro research use only