

Neurotrophin-3

Catalog No: #AP85911



Package Size: #AP85911-1 50ug #AP85911-2 100ug #AP85911-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neurotrophin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	NT-3
Accession No.	Uniprot ID:Q06AV0
Uniprot	Q06AV0
GeneID	100144493;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NTF3

Note: This product is for in vitro research use only