

Myosin-binding protein C, cardiac-type

Catalog No: #AP85913



Package Size: #AP85913-1 50ug #AP85913-2 100ug #AP85913-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myosin-binding protein C, cardiac-type
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cardiac MyBP-C,C-protein, cardiac muscle isoform
Accession No.	Uniprot ID:P56741
Uniprot	P56741
GeneID	295929;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mybpc3

Note: This product is for in vitro research use only