

Nuclear receptor subfamily 2 group C member 2

Catalog No: #AP85920



Package Size: #AP85920-1 50ug #AP85920-2 100ug #AP85920-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nuclear receptor subfamily 2 group C member 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Orphan nuclear receptor TAK1,Orphan nuclear receptor TR4,Testicular receptor 4,Mtr2r1,Tak1,Tr4
Accession No.	Uniprot ID:P49117
Uniprot	P49117
GeneID	22026;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nr2c2

Note: This product is for in vitro research use only