

Nuclear receptor subfamily 5 group A member 2

Catalog No: #AP85921



Package Size: #AP85921-1 50ug #AP85921-2 100ug #AP85921-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Nuclear receptor subfamily 5 group A member 2 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | Liver receptor homolog 1,LRH-1,Lrh1 |
| Accession No. | Uniprot ID:P45448 |
| Uniprot | P45448 |
| GeneID | 26424; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nr5a2

Note: This product is for in vitro research use only