

Nitric oxide synthase, brain

Catalog No: #AP85931



Package Size: #AP85931-1 50ug #AP85931-2 100ug #AP85931-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nitric oxide synthase, brain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Constitutive NOS,NC-NOS,NOS type I,Neuronal NOS,N-NOS,Peptidyl-cysteine S-nitrosylase NOS1,bNOS
Accession No.	Uniprot ID:P29475
Uniprot	P29475
GeneID	4842;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NOS1

Note: This product is for in vitro research use only