

Epididymal secretory protein E1

Catalog No: #AP85932



Package Size: #AP85932-1 50ug #AP85932-2 100ug #AP85932-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Epididymal secretory protein E1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Niemann Pick type C2 protein homolog,cE1
Accession No.	Uniprot ID:Q28895
Uniprot	Q28895
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name: NPC2

Note: This product is for in vitro research use only