

NEDD8 ultimate buster 1

Catalog No: #AP85935



Package Size: #AP85935-1 50ug #AP85935-2 100ug #AP85935-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	NEDD8 ultimate buster 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Negative regulator of ubiquitin-like proteins 1
Accession No.	Uniprot ID:Q8MJ87
Uniprot	Q8MJ87
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NUB1

Note: This product is for in vitro research use only