

Myomesin-3

Catalog No: #AP85939



Package Size: #AP85939-1 50ug #AP85939-2 100ug #AP85939-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Myomesin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Myomesin family member 3
Accession No.	Uniprot ID:Q5VTT5
Uniprot	Q5VTT5
GeneID	127294;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MYOM3

Note: This product is for in vitro research use only