

# Matrix metalloproteinase-14

Catalog No: #AP85940

Package Size: #AP85940-1 50ug #AP85940-2 100ug #AP85940-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Matrix metalloproteinase-14
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MMP-14,MMP-X1,Membrane-type matrix metalloproteinase 1,MT-MMP 1,Membrane-type-1 matrix metalloproteinase,MT1-MMP
Accession No.	Uniprot ID:P50281
Uniprot	P50281
GeneID	4323;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MMP14

Note: This product is for in vitro research use only