

Nuclear apoptosis-inducing factor 1

Catalog No: #AP85943



Package Size: #AP85943-1 50ug #AP85943-2 100ug #AP85943-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Nuclear apoptosis-inducing factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Accession No.	Uniprot ID:Q6PFD7
Uniprot	Q6PFD7
GeneID	71254;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Naif1

Note: This product is for in vitro research use only