

Muellerian-inhibiting factor

Catalog No: #AP85946



Package Size: #AP85946-1 50ug #AP85946-2 100ug #AP85946-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Muellerian-inhibiting factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Anti-Muellerian hormone,AMH,Muellerian-inhibiting substance,MIS
Accession No.	Uniprot ID:P03972
Uniprot	P03972
GeneID	280718;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AMH

Note: This product is for in vitro research use only