

PWWP domain-containing protein MUM1

Catalog No: #AP85951



Package Size: #AP85951-1 50ug #AP85951-2 100ug #AP85951-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PWWP domain-containing protein MUM1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Mutated melanoma-associated antigen 1,MUM-1
Accession No.	Uniprot ID:Q08DK9
Uniprot	Q08DK9
GeneID	513471;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MUM1

Note: This product is for in vitro research use only