

Arachidonate 15-lipoxygenase B

Catalog No: #AP85956



Package Size: #AP85956-1 50ug #AP85956-2 100ug #AP85956-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Arachidonate 15-lipoxygenase B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	15-LOX-B,15-lipoxygenase 2,15-LOX-2,Arachidonate 15-lipoxygenase type II,Linoleate 13-lipoxygenase 15-LOb
Accession No.	Uniprot ID:Q8K4F2
Uniprot	Q8K4F2
GeneID	266604;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Alox15b

Note: This product is for in vitro research use only