

Mitogen-activated protein kinase kinase kinase 5

Catalog No: #AP85962



Package Size: #AP85962-1 50ug #AP85962-2 100ug #AP85962-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|-----------------------------------------------------------------------------------------------------|
| Product Name | Mitogen-activated protein kinase kinase kinase 5 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | Apoptosis signal-regulating kinase 1,ASK-1,MAPK/ERK kinase kinase 5,MEK kinase 5,ASK1,MAPKKK5,MEKK5 |
| Accession No. | Uniprot ID:Q99683 |
| Uniprot | Q99683 |
| GeneID | 4217; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MAP3K5

Note: This product is for in vitro research use only