

Matrix metalloproteinase-17

Catalog No: #AP85963



Package Size: #AP85963-1 50ug #AP85963-2 100ug #AP85963-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Matrix metalloproteinase-17
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MMP-17,Membrane-type matrix metalloproteinase 4,MT-MMP 4,Membrane-type-4 matrix metalloproteinase,MT4-MMP,MT4MMP
Accession No.	Uniprot ID:Q9ULZ9
Uniprot	Q9ULZ9
GeneID	4326;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MMP17

Note: This product is for in vitro research use only