

Stromelysin-1

Catalog No: #AP85964



Package Size: #AP85964-1 50ug #AP85964-2 100ug #AP85964-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Stromelysin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SL-1,EMS-2,Matrix metalloproteinase-3,MMP-3,Transin-1
Accession No.	Uniprot ID:P28862
Uniprot	P28862
GeneID	17392;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mmp3

Note: This product is for in vitro research use only