

Neudesin

Catalog No: #AP85966



Package Size: #AP85966-1 50ug #AP85966-2 100ug #AP85966-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neudesin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Neuron-derived neurotrophic factor, Secreted protein of unknown function, SPUF protein, Spuf
Accession No.	Uniprot ID:Q9CQ45
Uniprot	Q9CQ45
GeneID	66208;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Nenf

Note: This product is for in vitro research use only