

NADH dehydrogenase [ubiquinone] iron-sulfur protein 7, mitochondrial

Catalog No: #AP85972

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP85972-1 50ug #AP85972-2 100ug #AP85972-3 1mg

Description

Product Name	NADH dehydrogenase [ubiquinone] iron-sulfur protein 7, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Complex I-20kD, CI-20kD, NADH-ubiquinone oxidoreductase 20 kDa subunit, PSST subunit
Accession No.	Uniprot ID:O75251
Uniprot	O75251
GeneID	374291;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NDUFS7

Note: This product is for in vitro research use only