

Methylenetetrahydrofolate reductase

Catalog No: #AP85975



Package Size: #AP85975-1 50ug #AP85975-2 100ug #AP85975-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Methylenetetrahydrofolate reductase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Accession No.	Uniprot ID:Q5I598
Uniprot	Q5I598
GeneID	497032;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MTHFR

Note: This product is for in vitro research use only