

Mitogen-activated protein kinase kinase kinase 5

Catalog No: #AP85977



Package Size: #AP85977-1 50ug #AP85977-2 100ug #AP85977-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mitogen-activated protein kinase kinase kinase 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Apoptosis signal-regulating kinase 1,ASK-1,MAPK/ERK kinase kinase 5,MEK kinase 5,Ask1,Mekk5
Accession No.	Uniprot ID:O35099
Uniprot	O35099
GeneID	26408;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Map3k5

Note: This product is for in vitro research use only