

NFAT activation molecule 1

Catalog No: #AP85978



Package Size: #AP85978-1 50ug #AP85978-2 100ug #AP85978-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NFAT activation molecule 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Calcineurin/NFAT-activating ITAM-containing protein,NFAT-activating protein with ITAM motif 1,Cnaip
Accession No.	Uniprot ID:Q8R4V1
Uniprot	Q8R4V1
GeneID	74039;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nfam1

Note: This product is for in vitro research use only