

NADH-cytochrome b5 reductase 2

Catalog No: #AP85980



Package Size: #AP85980-1 50ug #AP85980-2 100ug #AP85980-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NADH-cytochrome b5 reductase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	b5R.2,RCJMB04_32e2
Accession No.	Uniprot ID:Q5ZHX7
Uniprot	Q5ZHX7
GeneID	408182;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CYB5R2

Note: This product is for in vitro research use only