

Oviduct-specific glycoprotein

Catalog No: #AP85981



Package Size: #AP85981-1 50ug #AP85981-2 100ug #AP85981-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Oviduct-specific glycoprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Estrogen-dependent oviduct protein,Oviductal glycoprotein,Oviductin,Chit5,Ogp
Accession No.	Uniprot ID:Q62010
Uniprot	Q62010
GeneID	12659;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ovgp1

Note: This product is for in vitro research use only