

Prolyl 4-hydroxylase subunit alpha-1

Catalog No: #AP85983



Package Size: #AP85983-1 50ug #AP85983-2 100ug #AP85983-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prolyl 4-hydroxylase subunit alpha-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	4-PH alpha-1,Procollagen-proline,2-oxoglutarate-4-dioxygenase subunit alpha-1
Accession No.	Uniprot ID:Q1RMU3
Uniprot	Q1RMU3
GeneID	518288;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:P4HA1

Note: This product is for in vitro research use only