

N-acyl-phosphatidylethanolamine-hydrolyzing phospholipase D

Catalog No: #AP85987

Package Size: #AP85987-1 50ug #AP85987-2 100ug #AP85987-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	N-acyl-phosphatidylethanolamine-hydrolyzing phospholipase D
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	N-acyl phosphatidylethanolamine phospholipase D
Accession No.	Uniprot ID:Q58CN9
Uniprot	Q58CN9
GeneID	541291;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NAPEPLD

Note: This product is for in vitro research use only