

Neurofilament light polypeptide

Catalog No: #AP85991



Package Size: #AP85991-1 50ug #AP85991-2 100ug #AP85991-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neurofilament light polypeptide
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	NF-L,68 kDa neurofilament protein, Micro glutamic acid-rich protein, Neurofilament triplet L protein
Accession No.	Uniprot ID:P02548
Uniprot	P02548
GeneID	281348;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NEFL

Note: This product is for in vitro research use only