

Neural cell adhesion molecule 2

Catalog No: #AP85993



Package Size: #AP85993-1 50ug #AP85993-2 100ug #AP85993-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neural cell adhesion molecule 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	N-CAM-2,NCAM21
Accession No.	Uniprot ID:O15394
Uniprot	O15394
GeneID	4685;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NCAM2

Note: This product is for in vitro research use only