

Nitric oxide synthase, inducible

Catalog No: #AP85994



Package Size: #AP85994-1 50ug #AP85994-2 100ug #AP85994-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nitric oxide synthase, inducible
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Inducible NO synthase, Inducible NOS, NOS type II, NOSII, Peptidyl-cysteine S-nitrosylase NOS2
Accession No.	Uniprot ID:Q27995
Uniprot	Q27995
GeneID	282876;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NOS2

Note: This product is for in vitro research use only