

Neurogenic locus notch homolog protein 1

Catalog No: #AP85996



Package Size: #AP85996-1 50ug #AP85996-2 100ug #AP85996-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neurogenic locus notch homolog protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Notch 1, Motch A, mT14, p300, Motch
Accession No.	Uniprot ID:Q01705
Uniprot	Q01705
GeneID	18128;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Notch1

Note: This product is for in vitro research use only