

Nuclear factor of activated T-cells 5

Catalog No: #AP85997



Package Size: #AP85997-1 50ug #AP85997-2 100ug #AP85997-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nuclear factor of activated T-cells 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	NF-AT5, T-cell transcription factor NFAT5, Tonicity-responsive enhancer-binding protein, TonE-binding protein, KIAA0827, TONEBP
Accession No.	Uniprot ID:O94916
Uniprot	O94916
GeneID	10725;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NFAT5

Note: This product is for in vitro research use only