

Nuclear factor NF-kappa-B p100 subunit

Catalog No: #AP85998



Package Size: #AP85998-1 50ug #AP85998-2 100ug #AP85998-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nuclear factor NF-kappa-B p100 subunit
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	DNA-binding factor KBF2,H2TF1,Lymphocyte translocation chromosome 10 protein,Nuclear factor of kappa light polypeptide gene enhancer in B-cells 2,Oncogene Lyt-10,Lyt10,LYT10
Accession No.	Uniprot ID:Q00653
Uniprot	Q00653
GeneID	4791;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NFKB2

Note: This product is for in vitro research use only