

NLR family CARD domain-containing protein 4

Catalog No: #AP85999



Package Size: #AP85999-1 50ug #AP85999-2 100ug #AP85999-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NLR family CARD domain-containing protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Caspase recruitment domain-containing protein 12,Ice protease-activating factor,Ipaf,Card12,Ipaf
Accession No.	Uniprot ID:Q3UP24
Uniprot	Q3UP24
GeneID	268973;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nlrc4

Note: This product is for in vitro research use only