

Beta-microseminoprotein

Catalog No: #AP86002



Package Size: #AP86002-1 50ug #AP86002-2 100ug #AP86002-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-microseminoprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Prostate secreted seminal plasma protein,Prostate secretory protein of 94 amino acids,PSP-94,Psp94
Accession No.	Uniprot ID:O08540
Uniprot	O08540
GeneID	17695;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Msmb

Note: This product is for in vitro research use only