

Nuclear receptor coactivator 1

Catalog No: #AP86005



Package Size: #AP86005-1 50ug #AP86005-2 100ug #AP86005-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Nuclear receptor coactivator 1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Pig |
| Other Names | NCoA-1, Steroid receptor coactivator 1, SRC-1, SRC1 |
| Accession No. | Uniprot ID:Q4PJW2 |
| Uniprot | Q4PJW2 |
| GeneID | 574068; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NCOA1

Note: This product is for in vitro research use only