

Nuclear receptor coactivator 3

Catalog No: #AP86006



Package Size: #AP86006-1 50ug #AP86006-2 100ug #AP86006-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nuclear receptor coactivator 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	NCoA-3, Amplified in breast cancer-1 protein homolog, AIB-1, CBP-interacting protein, p/CIP, Receptor-associated coactivator 3, RAC-3, Steroid receptor coactivator protein 3, SRC-3, Thyroid hormone receptor activator molecule 1, ACTR, Aib1, Pcip, Rac3, Tram1
Accession No.	Uniprot ID: O09000
Uniprot	O09000
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0) gene_full_name: Ncoa3

Note: This product is for in vitro research use only