

Mucin-15

Catalog No: #AP86008



Package Size: #AP86008-1 50ug #AP86008-2 100ug #AP86008-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mucin-15
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MUC-15,UNQ750/PRO1481
Accession No.	Uniprot ID:Q8N387
Uniprot	Q8N387
GeneID	143662;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MUC15

Note: This product is for in vitro research use only