

Cell surface glycoprotein MUC18

Catalog No: #AP86009



Package Size: #AP86009-1 50ug #AP86009-2 100ug #AP86009-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Cell surface glycoprotein MUC18 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | Cell surface glycoprotein P1H12,Melanoma cell adhesion molecule,Melanoma-associated antigen A32,Melanoma-associated antigen MUC18,S-endo 1 endothelial-associated antigen,CD146,MUC18 |
| Accession No. | Uniprot ID:P43121 |
| Uniprot | P43121 |
| GeneID | 4162; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MCAM

Note: This product is for in vitro research use only