

## Nitric oxide synthase, brain

Catalog No: #AP86013



Package Size: #AP86013-1 50ug #AP86013-2 100ug #AP86013-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Nitric oxide synthase, brain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Constitutive NOS,NC-NOS,NOS type I,Neuronal NOS,N-NOS,Peptidyl-cysteine S-nitrosylase NOS1,bNOS
Accession No.	Uniprot ID:Q9Z0J4
Uniprot	Q9Z0J4
GeneID	18125;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Nos1

Note: This product is for in vitro research use only