

Protein N-terminal asparagine amidohydrolase

Catalog No: #AP86014



Package Size: #AP86014-1 50ug #AP86014-2 100ug #AP86014-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein N-terminal asparagine amidohydrolase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Protein NH2-terminal asparagine amidohydrolase,PNAA,Protein NH2-terminal asparagine deamidase,PNAD
Accession No.	Uniprot ID:Q96AB6
Uniprot	Q96AB6
GeneID	123803;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NTAN1

Note: This product is for in vitro research use only