

High affinity nerve growth factor receptor

Catalog No: #AP86015



Package Size: #AP86015-1 50ug #AP86015-2 100ug #AP86015-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	High affinity nerve growth factor receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Neurotrophic tyrosine kinase receptor type 1,Slow nerve growth factor receptor,p140-TrkA,Trk-A,Trk,Trka
Accession No.	Uniprot ID:P35739
Uniprot	P35739
GeneID	59109;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ntrk1

Note: This product is for in vitro research use only