

Nucleobindin-2

Catalog No: #AP86016



Package Size: #AP86016-1 50ug #AP86016-2 100ug #AP86016-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nucleobindin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	DNA-binding protein NEFA,Gastric cancer antigen Zg4,Prepronesfatin,NEFA
Accession No.	Uniprot ID:P80303
Uniprot	P80303
GeneID	4925;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NUCB2

Note: This product is for in vitro research use only