

Nucleolar and spindle-associated protein 1

Catalog No: #AP86017



Package Size: #AP86017-1 50ug #AP86017-2 100ug #AP86017-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nucleolar and spindle-associated protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	NuSAP,ANKT,BM-037,PRO0310
Accession No.	Uniprot ID:Q9BXS6
Uniprot	Q9BXS6
GeneID	51203;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NUSAP1

Note: This product is for in vitro research use only