

Propionyl-CoA carboxylase beta chain, mitochondrial

Catalog No: #AP86020



Package Size: #AP86020-1 50ug #AP86020-2 100ug #AP86020-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Propionyl-CoA carboxylase beta chain, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PCCase subunit beta, Propanoyl-CoA:carbon dioxide ligase subunit beta
Accession No.	Uniprot ID:Q2TBR0
Uniprot	Q2TBR0
GeneID	515902;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PCCB

Note: This product is for in vitro research use only