

Programmed cell death protein 4

Catalog No: #AP86021



Package Size: #AP86021-1 50ug #AP86021-2 100ug #AP86021-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Programmed cell death protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Protein MA-3,Topoisomerase-inhibitor suppressed protein,Ma3,Tis
Accession No.	Uniprot ID:Q61823
Uniprot	Q61823
GeneID	18569;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pdcd4

Note: This product is for in vitro research use only