

# cGMP-specific 3',5'-cyclic phosphodiesterase

Catalog No: #AP86022



Package Size: #AP86022-1 50ug #AP86022-2 100ug #AP86022-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	cGMP-specific 3',5'-cyclic phosphodiesterase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	cGMP-binding cGMP-specific phosphodiesterase,CGB-PDE,PDE5
Accession No.	Uniprot ID:O77746
Uniprot	O77746
GeneID	403825;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PDE5A

Note: This product is for in vitro research use only